



PRODUCT CATALOG

2021-2022
Ver.1



Artificial Teeth



What is AC?	2
NEW ACE ANTERIOR.....	3
NAPERCE POSTERIOR.....	3
EFUCERA AC POSTERIOR	3
MILLION POSTERIOR	4
FLAT AC POSTERIOR	4
What is FX?	6
FX ANTERIOR	7
FX POSTERIOR	7
What is NS?	9
CROWN NS ANTERIOR.....	10
EFUCERA NS POSTERIOR	10
What is PX?	12
CROWN PX ANTERIOR.....	13
SOLUUT PX ANTERIOR	13
EFUCERA PX POSTERIOR	14
COMBINATION SET Package	16
FULL SET Package	16
PCS Form Package	17
SHADE GUIDE AC	18
SHADE GUIDE NS	18
SHADE GUIDE PX	18
TEETH CABINET.....	18


What is AC?

Acrylic resins are widely used ingredients in artificial teeth manufacturing generated through chemical reaction by applying polymerization initiator and heat to a monomer. Derived from methyl methacrylate monomer (MMA), polymethyl methacrylate (PMMA) is a light material which does not significantly increase the weight of the denture and forms chemical bonding to a denture base as it is made of the same material. Having characteristic features of high translucency and ease of handling, it can be made into various shapes and shades.

In the modern era, people have become more health conscious and particular in choosing acrylic teeth suitable for their dental prosthesis, consequently we began to develop and supply high quality standard AC acrylic teeth products to meet the advancing market demand. Equipped with our decades of experience in artificial teeth manufacturing and very strict compliance with quality standards, we were able to meet these market demands. With primary focus on aesthetics, we have meticulously engineered each tooth's layer and gradation to successfully manifest the natural appearance in shape, shade and translucency. Our years of painstaking research and development and expertise in production process have led us to design the AC acrylic teeth with unparalleled resistance against everyday wear and tear.



On laboratory test trials, our AC acrylic teeth have shown outstanding resistance against stain and discoloration – thanks to our unique formulation and sophisticated polymerization technique which inhibits surface oxidation and tarnishing. All of these physical properties which are perfectly suited for functional dental prosthesis have given us the confidence to introduce our Japanese technology, AC acrylic teeth, onto the ever evolving dental market.



In pursuit of matching individual teeth shape, dimension and colour, we have customized a variety of moulds and shades readily available when required. Each of which are devised to naturally resemble and function like the real teeth. We have tailored to reproduce the physical essence of a smile by the combined aid of realistic mamelon and fluorescence effects. While teeth alignment is constructively harmonized to imitate the teeth-mouth feeling sensation, teeth occlusion is excellently corresponded to restore ideal mastication, improved chewing efficiency and enhanced denture stability. These attributes, together with its physical properties, have made our AC acrylic teeth recognized as the best choice in the dental market worldwide.

We hereby offer to you our competitive, well-known and globally trusted, high quality AC acrylic teeth.

NEW ACE ANTERIOR

Two-Layer Acrylic Resin Teeth



In full and partial denture cases, the resin teeth closely harmonize in shape and color with natural teeth and can be easily arranged, and the wax gum festooned without difficulty.

Upper 23 Moulds				
Basic Form	Mould			
Tapering	T1	T2	T3	T4
	T5	T6		
Tapering Long	T4	T5	T6	T7
Square	S2	S3	S4	S5
	S6	S7	S8	
Square Short	SS2	SS3		
Ovoid	O2	O3	O4	O5
Lower 12 Moulds				
Mould				
L2	L3	L4	L5	L6
L7	L8	L9	L10	L11
S3L	S4L			
Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5			
Packing	Upper Lower	6pcs / SET : 16SET / BOX		

NAPERCE POSTERIOR

Two-Layer Acrylic Resin Teeth



Mould	Upper / Lower					
	M28	M30	M32	M33	M34	M36
Shades	A1	A2	A3	A3.5		
	A4	B1	B2	B3		
	B4	C1	C2	C3		
	C4	D2	D3	D4		
	W0.5					
Packing	Upper Lower	8pcs / SET : 12SET / BOX				

The cusp angle of NAPERCE POSTERIOR is 30°.

EFUCERA AC POSTERIOR

Two-Layer Acrylic Resin Teeth



Mould	Upper / Lower				
	28	30	32	34	36
Shades	A1	A2	A3	A3.5	
	A4	B1	B2	B3	
	B4	C1	C2	C3	
	C4	D2	D3	D4	
	W0.5				
Packing	Upper Lower	8pcs / SET : 12SET / BOX			

The cusp angle of EFUCERA AC POSTERIOR is 20°.

MILLION POSTERIOR

One-Layer Acrylic Resin Teeth



Mould	Upper / Lower				
	28	29	30	31	32
Shades	A1	A2	A3	A3.5	
	A4	B1	B2	B3	
	B4	C1	C2	C3	
	C4	D2	D3	D4	
	W0.5				
Packing	Upper Lower	8pcs / SET : 12SET / BOX			

The cusp angle of **MILLION POSTERIOR** is 33°.

FLAT AC POSTERIOR

Two-Layer Acrylic Resin Teeth



Mould	Upper / Lower			
	30	32	34	
Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5			
Packing	Upper Lower	8pcs / SET : 12SET / BOX		

The cusp angle of **FLAT AC POSTERIOR** is 0°.


Combination Table					
NEW ACE ANTERIOR		NAPERCE POSTERIOR	EFUCERA AC POSTERIOR	MILLION POSTERIOR	FLAT AC POSTERIOR
Upper	Lower				
T1	L2	M30	28	29	-
T2	L2	M30	30	29	30
T3	L6	M32	30	31	30
T4	L4	M30	30	30	30
T5	L7	M32	32	31	32
T6	L7	M32	34	32	34
TL4	L6	M32 (M34)	34	31	34
TL5	L8	M33 (M34)	34	32	34
TL6	L9 (L8)	M33 (M34)	34	32	34
TL7	L11	M34	34	-	34
S2	S3L	M30	28	29	-
S3	S4L	M30	30	29	30
S4	L4	M32	32	30	32
S5	L5	M32	34	30	34
S6	L6	M32 (M34)	34	31	34
S7	L7	M34 (M36)	34	-	34
S8	L10	M36	36	-	-
SS2	S3L	M28	28	29	-
SS3	S4L	M30	30	29	30
O2	S3L	M28	28	29	-
O3	L3	M30	30	30	30
O4	S4L	M32	32	31	32
O5	L6	M32	34	32	34

What is FX?

It is widely known that conventional acrylic teeth are susceptible to abrasion. Acrylic teeth gradually wear down in the mouth over time. This process accelerates when the patient frequently eats abrasive foods. While maintaining the physical advantages of acrylic material, we made an attempt to improve the performance of our acrylic teeth by incorporating unique filler. Through extensive research, we identified all possible ingredients and variations, which were tested to withstand our manufacturing process and then subjected to laboratory trials. One filler demonstrated excellent performance and became part of our new formulation called FX.

Both AC and FX uses high quality acrylic material with the same degree of resistance property to stain and discoloration. The presence of fillers in acrylic polymer usually makes it susceptible to stains. However, as a result of the correct proportion of our unique filler, FX formulation negates the effects of staining agents.

We have also customized a variety of moulds and shades exclusively for the FX line that are readily available. These moulds, different in design to that of AC, NS and PX, offers a range of selection when a particular mould desired is cannot be found in AC, NS or PX teeth line.



Compared to conventional acrylic teeth, FX, with a hardness of Hv = 24, is stronger by as much as 20% against abrasion. This quality translates to stronger resistance against everyday wear and tear and therefore longer. FX is available in Eficera FX, 20 degree, and FX Posterior, 30 degree, to enhance chewing efficiency in a variety of cases.

We hereby offer to you our high performance, revolutionized FX acrylic resin teeth.



FX ANTERIOR

Two-Layer Highly Performed Acrylic Resin Teeth



Upper 16 Moulds				
Basic Form	Mould			
Tapering	T4	T5	T6	T7
Square	S4	S5	S6	S7
Square Short	SS4	SS5	SS6	SS7
Combination	C4	C5	C6	C7

Lower 8 Moulds			
Mould			
LA4	LA5	LA6	LA7
LB4	LB5	LB6	LB7

Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5			

Packing	Upper Lower	6pcs / SET : 16SET / BOX
---------	----------------	--------------------------

FX ANTERIOR is a full 3-D reproduction of natural teeth with improved labial ridge to emphasize the labial surface morphology. Arrangements duplicating natural teeth are possible.

FX POSTERIOR

Two-Layer Highly Performed Acrylic Resin Teeth



Mould	Upper / Lower					
	28	30	32	33	34	36
Shades	A1	A2	A3	A3.5		
	A4	B1	B2	B3		
	B4	C1	C2	C3		
	C4	D2	D3	D4		
	W0.5					

Packing	Upper Lower	8pcs / SET : 12SET / BOX
---------	----------------	--------------------------



The cusp angle of FX POSTERIOR is 30°.

Combination Table		
FX ANTERIOR		FX POSTERIOR
Upper	Lower	
T4	LB4	M28
T5	LB5	M30
T6	LB6	M32
T7	LB7	M33 (M34)
S4	LA4	M30
S5	LB6	M30
S6	LB7	M33 (M34)
S7	LA7	M33 (M34)
SS4	LB4	M28
SS5	LB5	M30
SS6	LA6	M32
SS7	LB7	M33 (M34)
C4	LA4	M30
C5	LA5	M30
C6	LA6	M33 (M34)
C7	LA7	M33 (M34)


What is NS?

For decades, we have been supplying the global dental market both with acrylic teeth and composite teeth. Throughout our experience, we have noticed that acrylic teeth users tend to seek acrylic teeth of higher quality than what they are using. While composite teeth users tend to seek alternative material of comparable quality, more affordable and resistant against staining agents. With this market need, we have searched for the most suitable material in order to fill the gap between conventional acrylic and composite teeth in terms of quality performance and competitiveness in the market.

Addressing the stain susceptibility issue of composite teeth, we have chosen to keep the acrylic nature of the desired artificial teeth material. While we also know that incorporating large amount of filler in the acrylic material to enhance its physical properties would make it susceptible to stains. Along with our years of research, we have found the right material of desired quality that has led us to the development of a new artificial teeth product line called hard acrylic NS.



Unlike AC or FX, embedded inside the NS are very minute particles called nanoSilica that made its polymer matrix structure more compact and tougher. These nano-sized Silica particles strengthen the bonding between polymer strands making it harder and resistant against abrasion. Possessing hardness of Hv = 25, performance test showed that NS is 60% stronger than conventional acrylic material against abrasion. Thus, NS has opened the opportunity for users, who are not quite satisfied with conventional acrylic resin teeth, a higher quality and competitive three-layer alternative choice.



The market demand for PX moulds at competitive level has been in our list for many years. This demand has made us to decide creating NS moulds the same as those of PX and made available in complete VITA shades.

We hereby offer to you new NS that will challenge the smile of the industry!

CROWN NS ANTERIOR

Three-Layer Hard Acrylic nanoSilica-Reinforced Resin Teeth



Upper 24 Moulds				
Basic Form	Mould			
Tapering	T41	T51	T61	
Tapering Short	T41S	T51S	T61S	
Square	S51	S71	S81	
Square Short	S41S	S42S	S43S	S44S
	S51S	S52S	S61S	
Ovoid	O41			
Ovoid Short	O31S	O51S	O61S	
Combination	C41	C42	C51	C61

Lower 8 Moulds				
Mould				
N31S	N61S	N31	N32	N41
N42	N81	N71L		

Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5			
Packing	Upper Lower	6pcs / SET : 16SET / BOX		

CROWN NS ANTERIOR is a 3D-digital reproduction of natural anterior teeth. It features solid moulds with supplementary labio-lingual width and emphasized tubercle protrusion to render space clearance provided for easy adjustments and strong clutching on the lingual gum, respectively.

EFUCERA NS POSTERIOR

Three-Layer Hard Acrylic nanoSilica-Reinforced Resin Teeth



Mould	Upper / Lower				
		28	30	32	34
Shades	A1	A2	A3	A3.5	
	A4	B1	B2	B3	
	B4	C1	C2	C3	
	C4	D2	D3	D4	
	W0.5				
Packing	Upper Lower	8pcs / SET : 12SET / BOX			




The cusp angle of EFUCERA NS POSTERIOR is 20°.

Combination Table		
CROWN NS ANTERIOR		EFUCERA NS POSTERIOR
Upper	Lower	
T41	N32	28
T51	N42	30
T61	N61S	34
T41S	N32	28
T51S	N42	30
T61S	N61S	34
S51	N42	30
S71	N71L	34
S81	N81	36
S43S	N41	28
S44S	N41	28
S41S	N32	28
S42S	N31	28
S52S	N42	30
S51S	N42	30
S61S	N61S	34
O41	N32	28
O31S	N31S	28
O51S	N61S	32
O61S	N61S	32
C41	N41	32
C42	N41	28
C51	N42	30
C61	N61S	34

What is PX?

Acrylic resin teeth are widely known for their beauty in shades and shapes despite of the fact that their surface property deteriorates through long time of use. In order to address this weakness, we have been supplying composite resin teeth in the dental market. While it is true that composite resin teeth are much harder than those of acrylics, which prove high endurance in clinical use, they are much susceptible to stains. Composite resin teeth, in general, consist of stain-causing components – Urethane dimethacrylate (UDMA) and/or Bisphenol A-glycidyl methacrylate (Bis-GMA) or Bis-GMA analog, and filler. Recognizing these inherent weaknesses of both acrylics and composites, we made an attempt to remediate this problem.

Our endeavor of producing high endurance and stain resistant resin teeth made-up of single composite material has been realized through the development of PX. Possessing a hardness of Hv = 45, PX is more than 5 times stronger against abrasion which translates in superior protection against wear and tear, and much longer life on usage compared to acrylic materials. Our PX is the hardest composite resin teeth around the world!



After testing stain-repelling agents that are compatible with our production process and PX formulation, one exceptional fluorine-containing monomer showed satisfactory results. This monomer acts as teeth surface shield against stain-causing agents, and thus protects the stain susceptible composite matrix. Through clinical testing it has been proven that PX is twice as hard as acrylics, while demonstrating a similar stain resistance capacity as acrylics. The superior qualities exhibited by PX guided us to advanced composite resin teeth technology.

We hereby offer to you the hardest and stain resistant composite resin teeth you have been looking for!



CROWN PX ANTERIOR

Three-Layer Composite Resin Teeth



CROWN PX ANTERIOR is a 3D-digital reproduction of natural anterior teeth. It features solid moulds with supplementary labio-lingual width and emphasized tubercle protrusion to render space clearance provided for easy adjustments and strong clutching on the lingual gum, respectively.

Upper 24 Moulds				
Basic Form	Mould			
Tapering	T41	T51	T61	
Tapering Short	T41S	T51S	T61S	
Square	S51	S71	S81	
Square Short	S41S	S42S	S43S	S44S
	S51S	S52S	S61S	
Ovoid	O41			
Ovoid Short	O31S	O51S	O61S	
Combination	C41	C42	C51	C61

Lower 8 Moulds				
Mould				
N31S	N61S	N31	N32	N41
N42	N81	N71L		

Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5			
Packing	Upper Lower	6pcs / SET : 16SET / BOX		

SOLUUT PX ANTERIOR

Three-Layer Composite Resin Teeth



With **SOLUUT PX ANTERIOR**, the cervical and incisal area of the Anterior are emphasized in order to render natural appearance and secured with sufficient dentin layers in order to avoid unnecessary translucency effect, respectively.

Upper 24 Moulds				
Basic Form	Mould			
Tapering	T4	T5	T6	T7
Square	S4	S5	S6	S7
Square Short	SS4	SS5	SS6	SS7
Ovoid	O4	O5	O6	O7
Combination	C4	C5	C6	C7
Combination SP	CSP4	CSP5	CSP6	CSP7

Lower 8 Moulds				
Mould				
L4	L5	L6	L7	
LS4	LS5	LS6	LS7	

Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5			
Packing	Upper Lower	6pcs / SET : 16SET / BOX		

EFUCERA PX POSTERIOR


Three-Layer Composite Resin Teeth



Mould	Upper / Lower				
	28	30	32	34	36
Shades	A1	A2	A3		A3.5
	A4	B1	B2	B3	
	B4	C1	C2	C3	
	C4	D2	D3	D4	
	W0.5				
Packing	Upper Lower	8pcs / SET : 12SET / BOX			


The cusp angle of EFUCERA PX POSTERIOR is 20°.

EFUCERA PX




Full Balanced Form

Occlusal ridges have been eccentrically positioned in order to achieve full occlusal equilibrium. Contact points, A, B and C, on the occlusal surface are designed for enhanced denture stability. Contact points, A and C, are reserved for lateral movements during mastication.



Setting Line

Imaginary setting line is reserved for the technician's own denture arrangement in response for patient's distinct requirement. Thus, this design is mostly applicable for partial dentures requiring unique arrangement with respect to its consequent natural teeth.



(Ideal Anatomic Proximate) IAP Face

To aesthetically harmonize the teeth arrangement from the anterior tooth to the molar part, the buccal side has one-tooth to two-tooth overlapping relationship.

Combination Table		
CROWN PX ANTERIOR		EFUCERA PX POSTERIOR
Upper	Lower	
T41	N32	28
T51	N42	30
T61	N61S	34
T41S	N32	28
T51S	N42	30
T61S	N61S	34
S51	N42	30
S71	N71L	34
S81	N81	36
S43S	N41	28
S44S	N41	28
S41S	N32	28
S42S	N31	28
S52S	N42	30
S51S	N42	30
S61S	N61S	34
O41	N32	28
O31S	N31S	28
O51S	N61S	32
O61S	N61S	32
C41	N41	32
C42	N41	28
C51	N42	30
C61	N61S	34

Combination Table		
SOLUUT PX ANTERIOR		EFUCERA PX POSTERIOR
Upper	Lower	
T4	L4	28
T5	L5	30
T6	L6	32
T7	L7	32
S4	L4	28
S5	L5	30
S6	L6	32
S7	L7	32
SS4	LS4	28
SS5	LS5	30
SS6	LS6	32
SS7	LS7	32
O4	LS4	28
O5	LS5	30
O6	LS6	32
O7	LS7	32
C4	L4	28
C5	L5	30
C6	L6	32
C7	L7	32
CSP4	L4	28
CSP5	L5	30
CSP6	L6	32
CSP7	L7	32

COMBINATION SET Package



Packing 28pcs / set : 4set / box

This package is available for all artificial teeth.

FULL SET Package



Packing 28pcs / set : 6set / box

This package is available for all artificial teeth.

PCS Form Package

What is Pieces Form?

New Bulk Package – Making Big Small. While not only pursuing improvements in the quality of our artificial teeth, we also focused on the most efficient for of packing to you give you more space and easy access. With Pieces Form, the teeth are now free from their plastic plate and can be picked out easily and quickly.


There are 6 cell for Anterior and 8 cells for Posterior and each cell contains 20 teeth.
All the information you need is indicated on the side label.



Teeth can be picked/shaken out through the opening in the lid. Turn the lid until the arrow points to the type you need. Then just shake out to dispense the tooth.



20 full conventional sets can now be stocked by piling 4 cases of Anterior Upper/Lower and Posterior Upper/Lower. This is more efficient way of stocking your teeth.



Once used, the containers can be refilled with our Refill-Pack offering a more economic, efficient and waste reducing system. The Refill-Pack contains 20 teeth per bag.



Teeth Formula

		UPPER															
		7	6	5	4	3	2	1	1	2	3	4	5	6	7		
RIGHT		7	6	5	4	3	2	1	1	2	3	4	5	6	7		LEFT
		LOWER															

	Case	Refill
Anterior	6 parts x 20 pcs each (120pcs/case)	1 part x 20pcs/pack
Posterior	8 parts x 20 pcs each (160pcs/case)	1 part x 20pcs/pack

SHADE GUIDE AC

Shade Guide for Acrylic Resin Teeth



SHADE GUIDE NS

Shade Guide for Hard Acrylic Resin Teeth



SHADE GUIDE PX

Shade Guide for Composite Resin Teeth



TEETH CABINET



Packing	1 Unit / 6-Pallet Drawer
Dimension	1 Unit (W285 x D310 x H220)mm

Each pallet has a capacity to accommodate 48 or 36 Yamahachi Anterior or Posterior sets, respectively.

CAD/CAM Milling Materials



ARTESANO	20
PMMA BLOCK (with Pin)	20
PMMA BLOCK (without Pin) for ROLAND DWX-4 ...	20
PMMA DISK	21
PMMA DISK ZZ	21
PMMA DISK AG	21
WAX DISK	22
WAX DISK α	22
WAX DISK ZZ	22
WAX DISK AG	23
WAX BLOCK (without Pin) for ROLAND DWX-4 ...	23

ARTESANO

CAD/CAM Milling Hybrid Composite Resin Block Material



Type	Block with pin				
Size	S (10 x 12 x 15 mm)		M (12 x 14 x 18 mm)		
Packing	5 pcs / box				
Shades	A1	A2	A3	A3.5	A4
Usage	Crown / Inlays				
Physical Properties					
3-Point Flexural Strength, MPa	195				
Biaxial Flexural Strength, MPa	230				
Compression Strength, MPa	526				
Vickers Hardness	71				
Fluorescence	Yes				

PMMA BLOCK (with Pin)


CAD/CAM Milling Acrylic Material



Type	Block with pin				
Size	S (15.4 x 19 x 39 mm)		M (15.5 x 19 x 55 mm)		
Packing	10 pcs / box				
Shades	A1	A2	A3	A3.5	
	A4	B1	B2	B3	
	B4	C1	C2	C3	
	C4	D2	D3	D4	
	W0.5	CLEAR			
	For shades other than A1, A2 and A3, minimum order quantity is 20boxes				
Usage	Temporary crowns and bridge				
	Model framework for casting				

PMMA BLOCK (without Pin) for ROLAND DWX-4

CAD/CAM Milling Acrylic Material



Type	Block for Roland DWX-4		
Size	76 x 40 x 20		
Packing	5 pcs / box		
Shades	A1	A2	A3
Usage	Temporary crowns and bridge		
	Model framework for casting		

PMMA DISK

CAD/CAM Milling Acrylic Material



Type	Open			
Diameter(mm)	98.5			
Thickness(mm)	10	12	14	
	15	16	18	
	20	22	25	
Packing	1 pc / box			
Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5	Clear	*V-Pink	
	V-Pink is vein fibers-containing shade, available in 20, 22, 25 and 30mm. Clear is also available in 30mm.			
Usage	Temporary crowns and bridges			
	Model framework for casting			
	Shade V-Pink: Denture Production			

PMMA DISK ZZ

CAD/CAM Milling Acrylic Material



Type	Zirkonzahn			
Diameter(mm)	95			
Thickness(mm)	15	16		
	18	20		
	22	25		
Packing	1 pc / box			
Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5	Clear	*V-Pink	
	V-Pink is vein fibers-containing shade, available in 20, 22, 25 and 30mm.			
Usage	Temporary crowns and bridges			
	Model framework for casting			
	Shade V-Pink: Denture Production			

PMMA DISK AG

CAD/CAM Milling Acrylic Material



Type	Amann Girrbach			
Diameter(mm)	101			
Thickness(mm)	13	20		
Packing	1 pc / box			
Shades	A1	A2	A3	A3.5
	A4	B1	B2	B3
	B4	C1	C2	C3
	C4	D2	D3	D4
	W0.5	Clear	*V-Pink	
	V-Pink is vein fibers-containing shade, available in 20mm.			
Usage	Temporary crowns and bridges			
	Model framework for casting			
	Shade V-Pink: Denture Production			

WAX DISK

CAD/CAM Milling Wax Material



Type	Open		
Diameter(mm)	98.5		
Thickness(mm)	10	12	14
	15	16	18
	20	22	25
Packing	1 pc / box		
Color	GREEN	IVORY	
Usage	Model framework for casting		

WAX DISK α

CAD/CAM Milling Wax Material



Type	Open		
Diameter(mm)	98.5		
Thickness(mm)	10	12	14
	15	16	18
	20	22	25
Packing	1 pc / box		
Color	Gray		
Usage	*Model framework invested with Cristobalite materials for rapid heating in castings of gold and palladium alloys. *Model framework for casting		

WAX DISK ZZ

CAD/CAM Milling Wax Material



Type	Zirkonzahn	
Diameter(mm)	95	
Thickness(mm)	15	16
	18	20
	22	25
Packing	1 pc / box	
Color	Green	Ivory
Usage	Model framework for casting	

WAX DISK AG


CAD/CAM Milling Wax Material



Type	Amann Girschbach	
Diameter(mm)	101	
Thickness(mm)	13	20
Packing	1 pc / box	
Color	Green	Ivory
Usage	Model framework for casting	

WAX BLOCK (without Pin) for ROLAND DWX-4

CAD/CAM Milling Wax Material



Type	Block for Roland DWX-4	
Size	76 x 40 x 20	
Packing	5 pcs / box	
Color	Green	
Usage	Model framework for casting	

Synthetic Resin



BASIS.....	26
BASIS HI	26
BASIS TWIN CURE	27
BASIS FLOW II.....	28
BASING RESIN and BASING RESIN α	28
RE-FINE BRIGHT	29
PROVIFINE	30
PATTERN BRIGHT	31
PARTIAL BRIGHT	31
ORTHO BRIGHT	32
ORTHO BRIGHT COLOR	32
BASIS PC	33
BASIS ELAST	34
ACRY PELLET	34

BASIS

Acrylic Resin for Denture Base



BASIS is strong and exceptionally durable heat-curing acrylic resin for denture bases. It is comprised of various sized particles which reinforce denture solidity and enhance the structure. An adaptable and aesthetically pleasing denture is achievable without any air bubble formation or shrinkage.



Heat-Curing Method: Immerse the flask in a container of tap water. Apply heat gradually for about 30 minutes until boil. Let the resin completely cure for 30 – 40 minutes in boiling water. Cool the flask for about 30 minutes at room temperature. Recover denture after cooling completely.

Packing	Powder	Liquid	
Retail	1kg, 3kg, 10kg	500mL (Pipette x 1), 1L, 17L	
Shades (O-Pink, V-Pink, LFPink and LF α are vein fibers-containing shades.)			
Clear	O-Pink*	V-Pink*	LFPink*
LF α *	Light Pink	Pink	Dark Pink
Physical Properties			
Powder/Liquid Mixing Ratio, g:mL	100 : 43		
Flexural Strength, MPa	94		
Flexural Modulus, MPa	2391		
Vickers Hardness, Hv	22.9		
Sorption, $\mu\text{g}/\text{mm}^3$	24		
Solubility, $\mu\text{g}/\text{mm}^3$	0.4		

BASIS HI

Acrylic Resin for Denture Base



BASIS HI is an acrylic elastomere, high impact resistance heat-curing resin for denture bases. An adaptable and aesthetically pleasing denture of excellent temperature stability is achievable without any air bubble formation or shrinkage. High impact resistance – guarantees worry-free application and use for patients and dental professionals.

Heat-Curing Method: Immerse the flask in a container of tap water. Apply heat until boil. Let the resin completely cure for 30 – 40 minutes (Curing time starts when the water with the flask has started to boil). Cool the flask for about 30 minutes at room temperature. Recover denture after cooling completely.

Features:

- **High Impact Resistance.** **BASIS HI** is a mixture of acrylic and elastomeric polymers exhibiting both the advantages of typical plastic and rubbery materials. These combined characteristics resulted in **BASIS HI**'s superb durability.
- **Excellent Temperature Stability.** **BASIS HI** acrylic elastomeric formulation is designed to withstand thermal effect during denture production and usage. Denture integrity against deformation due to thermal effect is preserved.
- **Non-Creeping.** The right proportion of copolymer's cross-linking prevents the tendency of the denture from slow deformation inside the mouth's stress. Longevity of denture is guaranteed.
- **Outstanding Color Stability.** The problem over color tarnishing and fading is prevented by **BASIS HI**'s stable copolymers' cross-linking.
- **Allows More Sufficient Working Time.** An adaptable and aesthetically pleasing denture is achievable without any air bubble formation or shrinkage in a less working time.

Packing	Powder	Liquid (Basis)	
Retail	1kg, 3kg, 10kg	500mL (Pipette x 1), 1L, 17L	
Shades (All shades are vein fibers-containing shades.)			
O-Pink	V-Pink	LFPink	LF α
Physical Properties			
Parameter	Value		
Powder/Liquid Mixing Ratio, g:mL	100 : 43		
Flexural Strength, MPa	112.7		
Flexural Modulus, MPa	2400		
Vickers Hardness	19.5		
Sorption, $\mu\text{g}/\text{mm}^3$	24		
Solubility, $\mu\text{g}/\text{mm}^3$	1.2		
Residual Monomer, wt%	0.7		

* Please use with the BASIS Liquid.

BASIS TWIN CURE

Heat Shock and Microwave-Curing Resin for Denture base



SS FRP Flask for Microwave-Curing

BASIS TWIN CURE is a denture base resin material applicable for both Heat Shock and Microwave-Curing methods. An adaptable, aesthetically pleasing and void-free denture is achievable without any air bubble formation or shrinkage in a less working time.


Heat Shock-Curing Method: Immerse the flask in boiling water for 15 minutes. Cool the flask for about 30 minutes at room temperature. Recover denture after cooling completely.

Microwave-Curing Method: Put the flask* into the microwave machine at 500W and cure for 3 minutes. In case where metal wire (clasp, etc.) is used, invest plaster and put water (about 180 mL) on the side of flask and then apply the microwave. Recover denture after cooling completely. *Use SS FRP microwave-curing flask.


Using conventional denture base resin, formation of void spaces translate into denture porosity thus prone to fractures, cracks and deformations. **BASIS TWIN CURE** eliminates void spaces formation that causes denture mechanical failures.

Packing	Powder	Liquid	
Retail	1kg, 10kg	500mL (Pipette x 1), 17L	
Shades (All shades are vein fibers-containing shades.)			
O-Pink	V-Pink	LFPink	LF α
Physical Properties		Curing Method	
Curing Method		Heat Shock	Microwave
Parameter	Value		
Powder/Liquid Mixing Ratio, g:mL	100 : 40		
Flexural Strength, MPa	89.3	81.4	
Flexural Modulus, MPa	2,167	2,273	
Sorption, $\mu\text{g}/\text{mm}^3$	23	23	
Solubility, $\mu\text{g}/\text{mm}^3$	0.1	0.1	
Residual Monomer, wt%	0.2	0.2	

Basis Twin Cure



Conventional Acrylic Resin



BASIS FLOW II

Multipurpose Self-Curing Pourable Acrylic Resin



Packing	Powder	Liquid
1-1Set	650 g	500 mL
	Accessories: (Plastic Cup, Spatula, Measuring Spoon, Cylinder Cup, Pipette) x 1 each	
Retail	500g, 10kg	500mL (Pipette x 1), 4L
*Shades (LFPink and LF α are vein fibers-containing shades.)		
	Clear	LFPink*
	LF α*	LF α*
Physical Properties		
Powder/Liquid Mixing Ratio, g:mL	100 : 60	
Flexural Strength, MPa	90.0	
Elasticity, mm	14.6	
Hardness, Hv	15.6	

BASIS FLOW II is a multipurpose pourable cold-curing acrylic resin that allows for sufficient working time and shortens total processing time.



Pressure-Curing Method: Pressure-polymerize the resin for 30 – 60 minutes at 55°C and 0.2MPa in a pressure pod.

BASING RESIN and BASING RESIN α

Self-Curing Acrylic Resin for Custom Trays and Base Plates



Types		
Product Name	Liquid Type	*Hardening Time, min
Basing Resin	Normal	5
	Slow	7
Basing Resin α	Normal	5
	Slow	7
Powder/Liquid Mixing Ratio, g:mL		100 : 35

* Hardening time value using prescribed powder/liquid mixing ratio at 23°C. Hardening time at lower and higher room temperature will become longer and shorter, respectively.

Packing	Powder	Liquid
1-1Set	1 kg	500 mL
	Accessories: Pipette x 1	
	Basing Resin α 1-1Set is not available.	
Retail	1kg, 10kg	500mL (Pipette x 1), 17L Basing Resin α available in 500mL only.
Shades	Pink	Blue

BASING RESIN and **BASING RESIN α** are self-curing, non-adhesive resins for base plates and individual trays. Non-adhesiveness offers moulding by spatula or fingers possible. **BASING RESIN α** is specially formulated for firmer adherence and easy handling of wax on bases and trays.

RE-FINE BRIGHT

Fast Setting Self-Curing Resin



1-1 Set Packing



Powder and Liquid Retail Packing

Packing	Powder	Liquid
1-1Set	250 g	260 mL
	Accessories: (Silicon Cup, Paint brush (Thin and Thick), Cylinder Cup, Pipette) x 1 each	
Retail	250 g	260 mL (Pipette x 1)
*Shades (O-Pink, V-Pink, LFPink and LF α are vein fibers-containing shades.)		
	Clear	O-Pink*
	LFPink*	LF α*
	V-Pink*	V-Pink*
	A2	A3
	A3.5	A3.5

Usage

- Production of inlays, temporary dental crowns and bridges
- Denture repairs

Physical Properties	
Parameter	Value
Powder/Liquid Mixing Ratio, g:mL	1 : 0.5
*Hardening Time (23°C)	3m 30s
Working Time	1m
Flexural Strength, MPa	75
Flexural Modulus, MPa	1,517
Vickers Hardness, Hv	11.4
Sorption, μg/mm ³	16
Solubility, μg/mm ³	2.2
Residual Monomer, wt%	3.3

* Hardening time value using prescribed powder/liquid mixing ratio at 23°C. Hardening time at lower and higher room temperature will become longer and shorter, respectively.

RE-FINE BRIGHT is a self-curing resin with excellent anti-discoloration properties – conventional problem of tarnishing is avoided and transparency retained. Component particles are of various sizes – strengthening, bonding and enhancing other physical properties. Superior shaving and cutting is possible. Enhanced operation – as desired cutting is achieved by uninhibited revolutions of bars and points, and exceptional mixing ability of the powder and liquids results in accurate reproductions.

PROVIFINE

Fast Setting Self-Curing Resin



PROVIFINE is a self-curing resin with improved physical properties.

Packing		Type	Physical Properties								
Powder	Liquid		Powder/Liquid Mixing Ratio, g:mL	**Hardening Time (23°C)	Working Time	Flexural Strength, MPa	Flexural Modulus, MPa	Vickers Hardness, Hv	Sorption, $\mu\text{g}/\text{mm}^3$	Solubility, $\mu\text{g}/\text{mm}^3$	Residual Monomer, wt%
50 g, 250 g	100 mL, 260 mL	Normal	100 : 50	4m 30s	1m 30s	82	1,750	14.8	16	2.3	1.7
		Fast		3m 30s	1m	89	1,920	15.1	16	2.3	1.8
Shades		Clear	LFPink	LF α	A1	A2	A3				

* (LFPink and LF α are vein fibers-containing shades.)
 ** Hardening time value using prescribed powder/liquid mixing ratio at 23°C. Hardening time at lower and higher room temperature will become longer and shorter, respectively.

Usage

- Production of inlays, dental crowns and bridges
- Denture repairs

Self-curing resin with High Liquidity at the Time of Pouring, Low Sagging, Easy to Build-Up at the Time of Brush Loading!
 Good Operability • Aesthetics • Durability, Suitable for Provisional Restoration

Operability

High Liquidity at the Time of Pouring



Due to high liquidity, the resin can be poured in to the fine details of the Silicone Core.

Low Sagging, Easy to Build-Up at the Time of Brush Loading



The resin is easy to shape without hanging and can be manipulated to the desired form.


In case where fluidity is desired and want to adapt the resin



After applying the powder and liquid on the brush tip, the resin will build up faster than usual.


Operation time can be selected according to liquid type to be used

Two types of liquid, normal and fast, with different curing times are made available. Preferred operation time can be selected. * Powder is common to normal and fast types liquid. Normal Type - Suitable when time clearance is required during work. Fast Type - Suitable for shortening the working time.



Aesthetics

Example Using Shade A3 for Repair



The color difference between pouring and brush loading methods is small, reproduction of stable color tone is possible. The color unevenness using brush loading method is reduced.

PATTERN BRIGHT

Self-Curing Acrylic Resin for Patterns



PATTERN BRIGHT is a self-curing resin for various pattern applications. With its very low polymerization shrinkage, as minimum as 0.72%, a compatible and satisfactory pattern is achieved. Hardening time is designed for speedy-work completion. When brush method is used, pattern production is made easy thanks to its excellent viscosity property. An almost no incineration residue results to smooth surface of the casting body, thus requires only minimal polishing.

Packing	Powder	Liquid	Shade
1-1Set	100 g	100 mL	Pink
	Accessories: (Silicon Cup, Paint Brush, Pipette) x 1 each		
Retail	100 g	100 mL (Pipette x 1)	
Usage	<ul style="list-style-type: none"> • Making out patterns on metal plates, lingual bars, palatal bars and connectors • Production of patterns on various clasps • Production of various Konuskronen telescope exterior crown patterns • Production of bonded bridge patterns • Temporary bonding of worn wax 		

Physical Properties		
Parameter	Value	
Powder/Liquid Mixing Ratio, g:mL	100 : 50	
*Hardening Time (23°C)	3m 20s	
Vickers Hardness, Hv	13.4	
Flexural Strength, MPa	60	
Post Setting Shrinkage, (23°C)	After 30 minutes	0.05%
	After 24 hours	0.06%
Polymerization Shrinkage (23°C)	0.72%	
Incineration Residue (700°C)	0.067%	

* Hardening time value using prescribed powder/liquid mixing ratio at 23°C. Hardening time at lower and higher room temperature will become longer and shorter, respectively.

PARTIAL BRIGHT

Self-Curing Acrylic Resin



1-1 Set Packing

Packing	Powder	Liquid	Shade
1-1Set	250 g	260 mL	*LFPink
	*LFPink is a vein fibers-containing shade. Accessories: (Silicon Cup, Paint Brushes (thin and thick), Pipette, Measuring Cup) x 1 each		
Retail	250 g	260 mL (Pipette x 1)	
Usage	<ul style="list-style-type: none"> • Partial Dentures • Copy Dentures • Denture Component Repairs • Dental Attachments 		

Physical Properties	
Parameter	Value
Flexural Strength, MPa	88.5
Dissolution, %	0.85
Stain (Fuchsin)	3.0

* Conditions: (Building-up Technique) Temperature 50°C water, Pressure: 2 atm, Polymerization Time: 30 minutes

Features:

- Sets and cures in about 12 minutes, allowing for ample time to mould.
- No air bubbles formed when applied, (when using either the building-up, pouring or spraying technique) making it very easy to handle.
- Silicone core and plaster core are included for measuring convenience. No need for a flask.
- Easily polished after setting, without burs and points being obstructed.



Powder and Liquid Retail Packing

ORTHO BRIGHT

Self-Curing Resin for Orthodontic Applications



Packing	Powder	Liquid
1-1 Set (Starter Kit)	100 g	70 mL
	Accessories: (Silicon Cup, Cylinder Cup, Powder Container, Pipette) X 1 Each; Pipette Nozzle X 3	
Retail	500 g	250 mL (Pipette x 1)
Shades		
	Clear	*Pink
	*Pink: The liquid is Pink.	

ORTHO BRIGHT COLOR

Self-Curing Resin for Orthodontic Applications



Packing	Powder	Liquid			
1-1 Set (Starter Kit)	50 g (x 5 shades)	70 mL			
	Accessories: (Silicon Cup, Cylinder Cup, Pipette, Shade Guide) X 1 Each; Pipette Nozzle X 3				
Retail	250 g	250 mL (Pipette x 1)			
Shades					
	Clear	Blue	Red	Orange	Green

ORTHO BRIGHT and ORTHO BRIGHT COLOR

Physical Property

Hardening Time and Hardening Process	Flow of Mixture in Wet Condition	Vicker's Hardness
Mud or sand-like state	17.0	13.5
Dough state	15.0	13.0
Rubbery, Cured state	13.0	13.3
ORTHO BRIGHT COLOR	11.0	12.5
T30	9.0	12.0
730	7.0	11.5
800	10.3	13.3
	ORTHO BRIGHT COLOR	ORTHO BRIGHT COLOR
	(cm)	(Hv)

Usage: All types of Splint, Functional Orthodontic Appliances, Deciduous Dentures, Temporary Dentures, Individual Trays

Features:

- Hardening time for complete polymerization reaction extends to about 8 minutes allowing for sufficient working time.
- Liquid monomer diffuses into the interstices of the polymer beads – releasing tension – migrates evenly and then absorbs by the matrix to form a homogenous fluid state. Diffusion of the liquid is like percolation of water into the sand. Excellent viscosity prevents the mixture fluid from sagging or slopping allowing for accurate control and shaping.
- Superior hardness ideal for orthodontic applications.

Methods of Use

Sprinkle Technique. Apply a separating agent for denture base to a plaster model. Perform preparation such as wax relief and fixing wires. Sprinkle liquid onto the powder until basement is formed. When the shine of the resin has disappeared, form the model using fingers. When resin elasticity is felt, immerse in water at 40-50°C (Placing in a pressure pot is recommended in order to minimize air bubble formation.)


Mixing Technique. Measure appropriate amount of powder and liquid. Put powder into liquid and mix using spatula or mixing stick. Mix slowly to avoid air bubble formation. When the mix has turned into paste-like body, pour into model. When the shine of the resin has disappeared, form the model using fingers. Use Sprinkle Technique for narrow parts. When resin elasticity is felt, immerse in water at 40-50°C (Placing in a pressure pot is recommended in order to minimize air bubble formation.)

Resin Packing Technique. Follow Mixing Technique for preparation. When the resin reaches the doughy state, immediately pack into the flask. Press the flask by hydraulic press until polymerization is complete (operate pressing before the curing process starts, refer to hardening time).

Brush On Technique. Put appropriate amount of powder and liquid to their corresponding containers. Wet the tip of the brush and dip into the powder. Take desired amount of powder to suffice powder load. Stack the load mixture until desired amount is achieved. Let hard-polymerize. Bigger brush is recommended for efficient results.

BASIS PC

Thermoplastic Resin Material for Denture Base (Polycarbonate)



Basis PC Pellets

Aluminum Tubes


BASIS PC is a new semi-flexible thermoplastic injection resin base material.

BASIS PC is allergic reaction-free, odorless and easy to polish. It is applicable for both full and partial dentures injection technique.

Features:


Balanced Strength

Basis PC mediates the gap between Nylon and Acrylic's strength characteristics resulting in its exceptional and distinctive quality.




Outstanding Stain Resistance

Basis PC effectively resists stain similar to that of acrylics. High ΔE^* (Reflectance) and ΔE^* (Transmittance) of Nylon indicates ineffective surface stain resistance and penetration of stain into the material, respectively (*The larger the ΔE , the more susceptible material is to stains.)



Superb Durability

Basis PC's hardness is proximate to that level of Nylon. Lower hardness, compared to acrylics, proves higher tenacity.




Legend:

- 1 Nylon
- 2 Basis PC
- 3 High Impact Acrylic
- 4 Conventional Acrylic


Semi-Flexible

The concurrent proportion of **Basis PC's** strength characteristics gave rise to a new semi-flexible denture base material.



Exceptional Water Sorption

Highly hydrophobic functionality resists water sorption. Propagation of bacteria-causing odor is controlled. Minimal deformation in intra-oral environment for long period of use is guaranteed.



Repairable


Material	Repair-ability*
Nylon	No
Basis PC	Yes
High Impact Acrylic	Yes
Conventional Acrylic	Yes

* Repair-ability using self-curing acrylic resin

BASIS ELAST

Thermoplastic Resin Material for Denture Base (Nylon)



Packing	Shade
300g, 1kg	Marble α 
Accessory	Aluminum Tubes (φ=2.5cm)
Physical Properties	
Parameter	Value
Flexural Strength MPa	83
Flexural Modulus, MPa	1392
Vickers Hardness, Hv	11.0
Elasticity, mm	8.7
Sorption, μg/mm ³	24.8


BASIS ELAST is a rigid-type and monomer-free (polyamido) Nylon denture base material with moderate elasticity suitable for non-metal clasp denture applications. **BASIS ELAST** is a flexible material with sufficient hardness for easy polishing.

Injection Parameters:
 Dry **Basis Elast** pellets at 80-90°C for 6hours before use;
 Melting Temperature 290°C; Melting Time 17min.; Injection Pressure 0.8MPa; Flask Temperature 60-90°C

ACRY PELLET

Thermoplastic Resin Material for Denture Base (Acrylic)



Packing	Shade
1kg	Marble H ※Vascular Pattern (without fiber) 
Physical Properties	
Parameter	Value
Flexural Strength MPa	75
Flexural Modulus, MPa	1810
Vickers Hardness, Hv	18
Sorption, μg/mm ³	20

High impact resistance and excellent toughness. The acrylic resin composition allows it to be used for repairing with self-curing resin and rebasing with relining materials.

Injection Conditions:
 Please set automatic oven at 80°C and use pellet after 6 hours of drying.
 Melting Temperature 275°C
 Dissolution Time 22 minutes
 Working Pressure 9atm
 Flask Heating 100°C

Waxes



PARAFFIN WAX	36
ROLLING WAX	36
CARVING WAX	36
BITE RIM STICK	36
DIPPING WAX	37
PRO UTILITY WAX	37
KOLBEN WAX	37
BITE WAX PRE-CUT TYPE	37
PRO LINE WAX	38

PARAFFIN WAX

Dental Use Paraffin Wax



Packing	Type	Color	Size
500g / 1kg / 5lbs	Soft / Medium	Light Pink / Pink	Regular (146 x 74 x 1.4mm) Large (170 x 85 x 1.4mm)

- Features:**
- Moderate plasticity and toughness.
 - Good crimping and retention of artificial teeth.

ROLLING WAX

Dental Use Sprue Wax



Packing	Diameter, mm	Shade
All sizes 270g / box	5	Pink
	6	
	7	
	8	

- Features:**
- User-friendly pliability for straightforward sprue applications
 - Convenient dispenser box guaranteeing wax protection

CARVING WAX

Dental Use Modeling / Waxing – up




Packing	Type		Shades		
	Cylinder	Stick	Ivory	Gray	
50g	140g (60 sticks)		Red	Blue	Green

- Features:**
- **Superb solidity.** Unaffected by varying atmospheric conditions
 - **High opacity and excellent color stability**
 - **Exceptional thermal expansion capacity.** Non-vulnerable to deformation due heat effects, robust shape guaranteed
 - **Burns out clean** with very little residue.
 - Minimal chipping, non-sticky to hands and instruments, outstanding shaving

BITE RIM STICK

Dental Use Pre-fabricated Wax for Occlusion Rims



Packing	Sizes*		Length	Shade
	S (Short)	L (Long)	25 cm	
	50 sticks / box (All sizes)			Pink

*sizes pertain to the arc length of the concavity

- Features:**
- Available in two sizes to appropriately fit the alveolar ridge's surface area
 - **No waste.** One stick sufficient for ridges of two full dentures

DIPPING WAX

Dental Use Coping Wax



Packing	Net Weight	Color
	200g	Yellow

Melting Range: (65 – 75) °C

- Features:**
- Optimum Elasticity
 - Burns out clean
 - Minimal Shrinkage
 - Excellent color stability even after repeat use

Relationship between Coping Thickness and Temperature

Temperature, °C	80	85	90
Thickness, mm	0.57	0.49	0.45

*Condition: Dipping Time 0.5second at 25°C

PRO UTILITY WAX

Dental Use Utility Wax



Packing	Sizes		Type / Hardness		Color
	Long (5x280)mm	Short (5x140)mm	Soft	Hard	
	125g / box				Red

- Features:**
- Soft, adhering and expandable wax
 - Soft and Hard types provide extensive range of practical applications
 - Ultimate variety in utility waxes

KOLBEN WAX

Dental Use Base Margin Forming Line Wax



Packing	Size	Color
	(2.2diameter x 200)mm	Red
	500pcs / box	

- Features:**
- Time-saving base margin and shape moulding wax
 - Easy to use and fix own design

BITE WAX PRE-CUT TYPE

Dental Use Pre-Cut Sheet Wax



Size	(137 x 73) mm
Pre-Cut Sheet Size	(15 x 73) mm
Packing	500g / box

- Features:**
- Wax for occlusion adjustment of natural teeth or denture.
 - Can be easily separated as they are pre-cut at 15mm- intervals.
 - Uses hard wax, minimal deformation can be achieved after bite-taking procedure.
 - Softens at low temperature, difficult to break even in the thin film state, can easily take the occlusion impression.
 - Occlusion impression is relatively easy to obtain with minimal strain and deformation.

PRO LINE WAX

Dental Use Pre-fabricated Casting Line Wax



Features:

- Exceptionally recommended for casting alloys for bases, clasps and sprue lines.
- **Optimum Elasticity.** High endurance over breaking on curve applications
- Superior welding abilities and applicable for wide range of uses
- **Mediated Casting Flow.** Glossy and smooth surface allows casting metal to flow easily

PRO LINE WAX Form and Packing					
Type	Shape	Diameter, mm	Height, mm	Usage	Packing, pcs / box
YR 05	•	(0.5)	-	Resin retaining Line of Metal Bases and Vents	120
YR 07	•	(0.7)	-		
YR 10	•	(1.0)	-		
YR 12	•	(1.2)	-		
YR 15	•	(1.5)	-	Sprue Line of Crowns, Bridges and Inlays	120
YR 20	•	(2.0)	-		
YR 25	•	(2.5)	-		
YR 32	•	(3.2)	-	Sprue Lines of Metal Bases	30
YR 35	•	(3.5)	-		
YR 40	•	(4.0)	-		
YR 50	•	(5.0)	-		
YR 60	•	(6.0)	-		10
YH 14	◐	1.4	1.1	Clasps	120
YH 16	◐	1.6	1.1		
YH 18	◐	1.8	1.1		
YH 19	◐	1.9	1.0		
YH 22	◐	2.2	1.2		
YH 28	◐	2.8	1.1		
YP I	◑	4.0	1.0	Palatal Bars	60
YP II	◑	4.0	1.5		
YL I	◒	3.1	1.4	Lingual Bars	60
YL II	◒	3.5	2.0		

Separating Agent and Cleansing Agent



APOLLON SEP (Normal)	40
APOLLON SEP (Low Viscosity).....	40
SPRAY BOTTLE	40
BREAK	40
WAX PATTERN CLEANER	40
WAX PATTERN CLEANER AQUA	40
APOLLON VARNISH	40
BRUSH CLEANER	40
TK SILICONE CLEANER	41
DOWEL PINS CLEANER	41
CLEAN UP.....	41
TRAY CLEANER (Powder).....	41
TRAY WASH (Liquid)	41
PIPE CLEAN (Liquid)	41
POLISH CLEANER	41
MIRROR CLEANER	41
HAND CLEANER	41

Separating Agent and Cleansing Agent

APOLLON SEP (Normal) Denture Base Separating Agent



Apollon Sep is a separating agent for resin denture bases with sodium alginate solution as the main ingredient, effective on flasks and plaster separation tasks.

Packing	500 mL	2 L	18 L
----------------	--------	-----	------

APOLLON SEP (Low Viscosity) Denture Base Separating Agent



Apollon Sep Low Viscosity offers easier work application. Handling becomes easier when used with the Spray Bottle.

SPRAY BOTTLE APOLLON SEP Accessory



Packing	300 mL (Empty)
----------------	----------------

*Spray Bottle is applicable only for Low Viscosity type.

BREAK Plaster Dissolution Agent



Packing	1 L	2 L	18 L
----------------	-----	-----	------

Dissolves plaster and gypsum left attached to dentures and cast materials. Progress of dissolution can be judged by the changing of the liquid color.

WAX PATTERN CLEANER Wax Pattern Strewing Agent



Packing	180 mL (Spray Type)
----------------	---------------------

Wax Pattern Cleaner application before investing enables for smooth painting of the investment and prevents porosity and uneven surface on the casting materials.

WAX PATTERN CLEANER AQUA Wax Pattern Strewing Agent (Spray-Type)



Packing	20 g / bottle
----------------	---------------

Lubricates casting surface, prevents bubble generation and uneven surface of the casting material. It can also be used for dental resin patterns since it does not contain ethanol.

APOLLON VARNISH Wax Pattern Separating Agent



Packing	100 mL (Paint Type)
----------------	---------------------

Wax pattern separating agent consisting mainly of surfactant for easy separation of applied wax pattern onto the surface of dentures, plaster casts and metals.

BRUSH CLEANER Brush Cleaner for Self-Curing Resin



Packing	100 mL
----------------	--------

Indication for Use:

- Removal of residual self-cure resin adhered on the brush
- Removal of polisher rouge stained on a casting object
- Removal of instant glue on a dowel pin

Usage: Pour appropriate amount of the liquid in a rubber cup, glass bottle or duppen glass. Immerse tissue paper for 5 minutes, and then use the wet tissue to wipe off the resins.

Separating Agent and Cleansing Agent

TK SILICONE CLEANER Silicone Surface Lubricating Agent



Packing	180 mL Spray Type (LPG)
----------------	-------------------------

Indication for Use:

- Spray TK Silicone Cleaner for smooth flow of model agent and prevention of bubbles from entering into the silicone impression or duplicate impression.

DOWEL PINS CLEANER Instant Glue Powerful Solvent



Packing	300 mL
----------------	--------

Dip the Dowel Pin with adhered instant glue into the ultrasonic cleaner with undiluted Dowel Pin Cleaner for 4 to 5 minutes.

When contact with fingers or hands, rub for 3 to 4 minutes with infiltrating absorbent cotton and wash using cold water.

CLEAN UP Non-Heating Gold and Palladium Alloys Cleaning Liquid



Packing	500 mL
----------------	--------

Clean Up is a cleaning agent for the removal of Gold oxide and Palladium oxide layers without evolution of heat. Please use undiluted liquid.

TRAY CLEANER (Powder) Alginate Impression Materials Cleaning Agent



Packing	1 kg (spoon included)
----------------	-----------------------

Tray Cleaner is a fast-acting tray cleaner for the removal of alginate impression materials by carbonization and simultaneously sterilizes and deodorizes the tray.

Usage: Mix 50g-100g of powder and mix with 1L of water.

*The powder dissolves faster at higher temperatures.

TRAY WASH (Liquid) Alginate Impression Materials Only Cleaning Agent



Packing	1000 mL
----------------	---------

Tray Wash is for rapid removal of alginate impression material adhered to trays. It is an excellent corrosion resistance agent for aluminum, nickel and chromium-plated trays.

Usage: Dilute with water by 10 parts. For severe dirt application, please dilute with water by 5 parts.

PIPE CLEAN (Liquid) Dental Drain Pipes Cleaner



Packing	1000 mL
----------------	---------

Pipe Clean has an excellent sterilizing and deodorizing abilities, it prevents the outbreak of unpleasant odors. It assists in washing off and decomposition of organic residues (blood, saliva, etc.) which can stain drainpipes and cuspidors.

Usage: Dilute with water by 10 parts. For severe dirt application, please dilute with water by 5 parts.

POLISH CLEANER Cleaning Liquid Exclusively for Ultrasonic Cleaners



Packing	1000 mL
----------------	---------

Polish Cleaner is developed as a cleaning agent for ultrasonic cleaners. It is transparent, rapidly removes all adhered rouge abrasives on the prosthetic appliances and eugenol cements.

MIRROR CLEANER Rouge-Type Abrasives Cleaning Liquid



Packing	1000 mL
----------------	---------

Mirror Cleaner is cleaning agent for rouges adhered to prosthetic appliances.

Usage: Mix 5mL of Mirror Cleaner with hot water for resins; mix 100mL of Mirror Cleaner with lukewarm water for metals and use ultrasonic cleaner for 2 - 3 minutes.

HAND CLEANER Hand Wash Powder Soap



Packing	1 kg
----------------	------

Hand Cleaner has an outstanding effect for washing hands after polishing works. It thoroughly cleans the dirt, sand and abrasives; it can also be used for cleaning various instruments.

Plaster and Investment Materials



FINE STONE.....	44
FINE ROCK	44
DENTAL PLASTER	44
CRISTO HEAT SHOCK	44
NEO WHITE	44

Plaster and Investment Materials

FINE STONE

Plaster for Hard Models



- Features:**
- Soft impression material for excellent reproductions.
 - Lighter shade allows for extended working time – without eye-strain.

Technical Data (23±5°C)		Packing	Shade
Mixing Ratio (Powder / Liquid)	100g / 24mL		
Hardening Time	10 min		
Hardening Shrinkage	0.16%		
Compression Strength	40 MPa		

FINE ROCK

Plaster for Hard Models



- Features:**
- Insignificant expansion ratio during setting.
 - Strength and hardness above standard.
 - Suitable for pouring into fine details.

Technical Data (23±5°C)		Packing	Shade
Mixing Ratio (Powder/Liquid)	100g / 22 mL		
Hardening Time	9 - 14 min		
Hardening Shrinkage	0.06 - 0.09%		
Compression Strength	≥ 40 MPa		

DENTAL PLASTER

Plaster for Dental Use



- Features:**
- High degree of purity, strength and consistent high quality.
 - Easily moulded with excellent precision during production achievable.

Standard Dental Plaster Technical Data (23±5°C)		Packing	Dental Technician Plaster Technical Data (23±5°C)		Packing
Mixing Ratio (Powder/Liquid)	100g / 40mL		10kg x 2 (20kg)	Mixing Ratio (Powder / Liquid)	
Hardening Time	12 min	Hardening Time		12 min	
Hardening Shrinkage	0.28%	Hardening Shrinkage		0.28%	
Compression Strength	17 MPa	Compression Strength		17 MPa	

CRISTO HEAT SHOCK

Cristobalite Investment Materials for Rapid Heating



Accurate casting investment material for gold and palladium alloys.

Technical Data (23±5°C)	
Mixing Ratio (Powder / Water)	100 g / 35mL
Hardening Time	15 min
Hardening Expansion	0.5 %
Thermal Expansion	1.4 %
Compression Strength	4.0 MPa
Packing	3kg

NEO WHITE

Dental Phosphate-Bonded Casting Investment Material



Technical Data (23±5°C)	
Mixing Ratio (Powder/Liquid)	100 g / 15 mL
Hardening Time	11 min
Average Use time	4 min
Hardening Expansion	0.75 %
Thermal Expansion	1.53 %
Heating Compression Strength	9 MPa

Packing		
	Powder	5 kg
Liquid	500 mL	

Uses: Partial Denture Bases, Bars and Clasps

- Features:**
- Minimal burning on casting surface, post-casting polishing dramatically reduced.
 - Casting easily cut from investment, therefore no damage for casting.
 - Excellent performance, fluidity and result reproducibility are outstanding.
 - Sufficient expansion confirmed excellent adaptability.

Precious Alloys, Non-Precious Alloys, Alloy Wires



NEORIUM S (Soft)	46
NEORIUM H (Hard)	46
NEO TITAN WIRE	47
FINE COBALT CLASP WIRE	47
YAMAHACHI SEMI CIRCLE WIRE (Regular/Soft)	47
REINFORCEMENT WIRE	47
YAMAHACHI LINGUAL BAR WIRE	47
YAMAHACHI PALATAL BAR WIRE	47
YAMAHACHI CLASP WIRE	47
LINGUAL BAR WIRE (Thin Type)	47

NEORIUM S (Soft)

Dental Casting Cobalt Chrome Alloy
(Exclusively for High Frequency Casting Machines)



Packing	1 kg / bottle	
Type	Coin	5g / pc
	Cylinder	10g / pc

Technical Data JIS T 6115

Metal	Co	Cr	Mo	Si, Mn, C, N, B*
Composition, %	65.2-59.8	28.0-30.0	5.6-5.9	1.2-4.3

*Others

Physical Properties	Value
Liquidus (Melting) Point	1,394 °C
Solidus Point	1,360 °C
Tensile Strength, MPa	≥ 685
Elongation, %	≥ 3
Hardness, Hv	≥ 340

NEORIUM H (Hard)

Dental Casting Cobalt Chrome Alloy
(Exclusively for High Frequency Casting Machines)



Packing	1 kg / bottle	
Type	Coin	5g / pc
	Cylinder	10g / pc

Technical Data JIS T 6115

Metal	Co	Cr	Mo	Si, Mn, C, N, B*
Composition, %	64.3-59.2	28.0-30.0	6.7-7.1	1.0 ~ 3.7

*Others

Physical Properties	Value
Liquidus (Melting) Point	1,385 °C
Solidus Point	1,355 °C
Tensile Strength, MPa	≥ 685
Elongation, %	≥ 3
Hardness, Hv	≥ 340

NEORIUM S and NEORIUM H

Cautions:

- For use in Argon Gas Atmospheric Melting Chambers only
- Not for use in Arc Casting Chambers

Uses: Full Denture Bases, Partial Denture Bases, Bars and Clasps

Features:

- Difficult to break, flexible casting achievable. Therefore the amount of adjusting to prevent casting defects is greatly decreased.
- Extractability from the investment material is excellent. Especially effective when used with Yamahachi investment **Neowhite**, the casting is easily removed from the investment material.

- **Hard to break even if casting deformation is adjusted.**

Due to sufficient elongation property, production of supple and hard to break casting is possible.

- **Neorium is made from powder metal ingredients.**

Compared to the dissolution method of production, the powder sintering method of production improves the alloy's physical properties because it utilizes more nitrogen and contained stably.

- Less deterioration even after reuse.

Neorium S	Virgin Material	First Reuse	Second Reuse
Elongation, %	13.8	11.5	12.1
Hardness, Hv	360	363	362
Tensile Strength, MPa	902	855	863

Neorium H	Virgin Material	First Reuse	Second Reuse
Elongation, %	8.8	8.5	8.6
Hardness, Hv	401	397	399
Tensile Strength, MPa	928	879	907

- Pellets are available in coin and cylindrical shapes.

NEO TITAN WIRE

Titanium Alloy Wire for Dental Use



Features:

- Ideal for areas with deep undercuts
- Superb corrosion resistance
- Excellent yield strength


Packing	2m / Roll			
Diameter, mm	0.7	0.8	0.9	1.0

Technical Data	Metal	Ti	Mo	Sn
Composition, %	80.5	10.8	8.7	

NEO TITAN WIRE 3 Main Advantages

1. Flexibility. Having 15 times elasticity limit than normal, NEO TITAN WIRE has higher limit against breakage.

	Neo Titan	Fine Cobalt*
Elasticity (Spring Property)	3.0%	0.2%



2. Low Allergy Risk. NEO TITAN WIRE is composed of elements (Ti, Mo, Sn) that have low toxicity and allergy risk.

Cellular Toxicity Low Limits	Ti Mo Sn Zr Nb Ta Pt
Cellular Toxicity Low Limits	Ni V Fe Co
Allergy Risk	Hg Ni Al Cd Cr Cu

3. Low Specific Gravity. Light NEO TITAN alloy results in minimum effect on the remaining teeth.

	Neo Titan	Fine Cobalt*
Specific Gravity	5.2	8.5

FINE COBALT CLASP WIRE

Dental Cobalt-Chromium Alloy Wire



Features:

- Exceptional elasticity, viscosity and corrosion resistance – all of which are essential qualities for clasp wire.
- Soldering easily accomplished.

Packing	5m / Roll		
Diameter, mm	0.8	0.9	1.0

Technical Data	Metal	Composition, %	Element	Composition, %
	Co	≥ 40.0	Mo	5.8 – 6.8
	Cr	20.5 – 22.5	Mn	0.9 – 1.5
	Ni	15.5 – 17.5	Si	≤ 0.5
	Fe	residual	C	0.10 – 0.15

REINFORCEMENT WIRE

Dental Stainless Steel Wire



Packing	6 m / Roll		
Flat Type	Hard / Soft		

Sizes	Width, mm	Height, mm
Thick	2.0	0.7
Medium	1.8	0.6

Semi-Circular Type	Hard / Soft		
Sizes	Width, mm	Height, mm	
Thick	2.0	1.0	
Medium	1.8	0.9	
Thin	1.4	0.7	

YAMAHACHI CLASP WIRE


Dental Stainless Steel Wire



Shape	Size	Shape	Diameter φ mm	Packing
Circular	0.8	○	0.8	5m Roll
	0.9	○	0.9	
	1.0	○	1.0	
	1.1	○	1.1	
	1.2	○	1.2	
	1.3	○	1.3	

YAMAHACHI SEMI CIRCLE WIRE (Regular/Soft)

Dental Stainless Steel Wire



Packing	3m / Roll		
----------------	-----------	--	--

Semi-Circular Type	Sizes	Diameter, mm	Height, mm
	1.4 – 1.4S	1.4	0.7
	1.6 – 1.6S	1.6	0.8
	1.8 – 1.8S	1.8	0.9
	2.0 – 2.0S	2.0	1.0
	2.3 – 2.3S	2.3	1.2

YAMAHACHI LINGUAL BAR WIRE YAMAHACHI PALATAL BAR WIRE

Dental Stainless Steel Wire



Packing	3pcs	
Length, cm	31	

Yamahachi Lingual Bar Wire	Sizes	Shape	Width, mm	Thickness, mm
	Small	□	2.2	1.2
	Medium	□	2.5	1.5
	Large	□	2.7	1.4

Yamahachi Palatal Bar Wire	Sizes	Shape	Width, mm	Thickness, mm
	Small	⤴	2.8	1.3
	Medium	⤴	3.0	1.4
	Large	⤴	3.7	1.25

LINGUAL BAR WIRE (Thin Type)

Dental Stainless Steel Wire



Packing	1m / Roll		
Sizes	Shape	Width, mm	Height, mm
S		3.0	0.9
SS		2.5	

Abrasive Materials / Polishing Materials



BRAZING DIA HP.....	50
CFP HOLDER	50
CERAMIC FIBER POINT	50
SILICONE BIG	51
TWISTER WHEEL	51
NEW SILICONE POINTS II	52
MANDREL CYLINDERS	52
URETHANE BIG	52
URETHANE DISK	52
ART POLISHER	53
YAMAHACHI CUTTING DISK	53
DIAMOND BRUSH	54
HOG(High Quality) HAIR BRUSH	54
HOG HAIR BRUSH	54
MIRROR BUFF.....	55
MILLION BUFF	55
MANDRELS #303	55
MP POWDER	55
MP BUFF	55
CREAMY SAND	55
SULFONE SAND	55
GLASS BEADS.....	55
ALUMINOUS.....	56
POLISHING POWDER	56
GRAZE POWDER	56
SILKY SHINE	56
BLUE SHINE.....	56
TIGER MULTI	57
TIGER MULTI MINI	57
TIGER MULTI GOLD	57
ARTE SHINE.....	57

BRAZING DIA HP

Dental Use Diamond Polisher

Number	#1	#2	#3	#4
Shape				
Number	#5	#6	#7	#8
Shape				

Uses: Modification Polisher for Porcelain

Features:

- **Excellent Polishing Ability.** Polishing surface made with sharp diamond grains.
- **High Clogging Resistance.** Diamond grain and physical object have large surface contact.
- **High Durability.** High chemical stability and mechanical retention of diamond grains during brazing.

Packing	Type	Coarseness	Working Speed	Color Code
1pc	1, 2, 3, 6, 7	Medium, Fine	Max. 30,000 rpm	Blue / Medium
	4, 5, 8	Fine		Red / Fine

CFP HOLDER

Dental Use Mandrel



Features:

- The shortened ceramic fiber can be extended by mounting in CFP Holder.
- * Please use glue when mounting the point in the holder.

Packing
5 pcs / case

CERAMIC FIBER POINT

Dental Use Polisher



Type	Green (Coarse)	#140 - #180
	Orange (Medium)	#320 - #400
	Red (Fine)	#1200 - #2000
Dimension	φ 2.35 x 50mm	
Packing	1pc / pack	

Uses:

- Polishing around pit fissures of inlay crowns
- Removal of air bubble inside the crown or clasp
- Fine adjustment or modification of resin and metal base or attachment
- Shape modification of Porcelain

Features:

- Sharp alumina fiber always protrudes on the surface allowing for excellent abrasion.
- Alumina fiber filled in high density packing to achieve clogging and minimal heat emission.
- Uniform-sized Alumina fibers packed in high density for reduced consumption.
- Does not break even at thinner diameter because of balanced required elasticity.

Attention: Operate at less than 20,000rpm. Follow the instruction of the hand-piece machine and check if the material is properly fixed. Check if material revolves evenly before use. Wear eye protector, mask for safe use. Do not use the product other than indicated by the manual.

SILICONE BIG

Silicone Big Points



Features:

- Excellent durability and stability.
- Wear occurs evenly and slowly, polishing power is exceptional.
- Seven varieties are available
- Unrivalled cost performance.

Packing	10pcs / box	100pcs / box
Working Speed	Max. 15,000rpm	
Size	L x W = 23 x 9.5 mm	

Specifications							
Appearance							
Polishing Code	C - 2	C - 3	M - 1	M - 2	M - 3	F - 2	F - 3
Polishing Texture	Coarse		Medium			Fine	
Color	Black	Gray	Dark Brown	Brown	Light Brown	Green	Light Green
Usage	Lab use composite	Acrylic	Amalgam, Precious Alloys	Lab use composite	Acrylic	Lab use composite	Acrylic

TWISTER WHEEL

Silicone Wheels



Features:

- Excellent durability and stability
- Slow wearing with exceptional polishing results; excellent cost performance

Packing	20pcs / box	50pcs / boxes
Working Speed	Max. 20,000rpm	
Size	All codes: D x W = 22 x 3.2 mm	

Specifications				
Appearance	Polishing Code	Polishing Texture	Color	Usage
	C - 2	Coarse	Black	Amalgam, Precious Alloys, Acrylic
	C - 3		Gray	Precious Alloys, Porcelain
	M - 1	Medium	Dark Brown	Cobalt-Chromium, Non-Precious Hard Alloys
	M - 2		Brown	Amalgam, Precious Alloys, Acrylic
	M - 3		Light Brown	Precious Alloys, Porcelain
	F - 2	Fine	Green	Amalgam, Precious Alloys, Acrylic
	F - 3		Light Green	Precious Alloys, Porcelain

NEW SILICONE POINTS II

Silicone Polisher



Polisher Code	Working Speed	Packing
#10, #13, #13S #28, #114	Max. 30,000rpm	12pcs / box 72pcs / box
#162	Max. 15,000rpm	12pcs / box 72pcs / box
Cylinder (Mandril x 1pc)		72pcs / box
Cup	Max. 30,000rpm	12pcs / box, 72pcs / box
Color	Polishing Texture	
Brown	Medium	
Green	Fine	

Specifications								
Shape								
Polisher Code	#10	#13	#13S	#28	#114	#162	Cup	Cylinder

Uses: Intermediate Polisher for Metal Alloys, Palladium Alloys, Acrylic Resin
Note: Cylinder-type available in Brown only

Features:
 • Contains combination of fine abrasive grains for shiny polish.

ART POLISHER

Silicone Wheel for Cobalt-Chrome Modifications



Packing	20pcs / box	50pcs / box	Polishing Usage
Working Speed	Max. 20,000rpm		
Size	D x W = 22.0 x 3.2 mm		
Type / Color	Hard / Blue (Medium Polish)	Soft / Light Blue (Medium Polish)	Cobalt-Chrome Au, Ag, Pd Alloys

YAMAHACHI CUTTING DISK

Metal Alloys Sprue Cutting Disks



Type	Size (Diameter x Thickness)	Packing	Usage	Working Speed
A	25 x 0.35 mm	50 pcs / box	Metal Alloy Sprues	Max. 15,000 rpm
B	25 x 0.60 mm	25 pcs / box	Metal Alloys	
C	38 x 0.60 mm	25 pcs / box	Metal Alloys	
E	22 x 0.23 mm	50 pcs / box	Ceramic	

Made with sharp edge to speedily cut sprues of silver, palladium alloy of course, nickel chrome alloy, up to cobalt-chrome alloy.

URETHANE DISK

Urethane Wheels



Packing	20pcs / box
Working Speed	Max. 15,000rpm
Size	D x W = 22 x 3.2 mm
Color	Polishing Texture
Blue (#100)	Coarse
Green (#320)	Medium

MANDREL CYLINDERS


Dental Use Mandrels



Features:	Packing
• Mandrels for Hand Piece Use New Silicone Point II and Cylinder Type Polishers	12 pcs / pack

URETHANE BIG

Urethane Big Points



Packing	10pcs / box
Working Speed	Max. 15,000rpm
Size	L x W = 24.5 x 9.5 mm

Features:
 • **Wobble-Free Polish.** Stable rotation and fine cushion from advanced Japanese technology result in ultra-smooth polishing experience.
 • **Efficient Bubble Buffer.** Heat-absorbing sponge-like polisher allows for heat-guarded and extended wear polishing.
 • **Multi-Purpose Polisher.** Highly effective polisher for wide range of applications: soft lining materials, mouthguards, splints, nylon, acrylic resin and metals.

DIAMOND BRUSH

Coarse Polishing Brush for Acrylic and Sulfonamide Resin



Features:

- Fiber brush is made up of specially formulated chemical fiber material that is static electricity inert – does not become dusty during polishing.
- Brush contains polishing powder material for fine polishing performance.
- Highly durable.

Packing	3 pcs / box	
Type	Soft (φ 67x15mm) Brush Line = 2	Regular (φ 67x15mm) Brush Line = 2 & 3

HOG HAIR BRUSH

Dental Polisher Hog Hair Lathe Brush



Features:

- Finest quality hog hair used making it suitable for coarse polishing of acrylic resin.
- Very satisfactory polishing performance is achieved when used with Sulfone Sand.

Packing	12 pcs / box		
Type	Number of Brush Line		
	1	2	3

HOG(High Quality) HAIR BRUSH

Dental Polisher Horse Hair Lathe Brush



Features:

- Center hub is made of solid wood resulting in minimal bristle loss.
- Bristle is made up of fine elastic material to assure good polishing contact resulting in excellent cleaning.

Packing	12 pcs / box		
Type	Number of Brush Line		
	1	2	3

MIRROR BUFF

Dental Polisher Finishing Buff



Features:

- Made from natural hemp suitable for finish polishing of metals and resins. High polishing capability and economical.
- Can skip sand paper process to cut down work time 3-5 times more efficient.

Packing	1 pc / pack
Size (Diameter x Thickness)	75 x 10 mm

MILLION BUFF

Dental Polisher Finishing Buff



Features:

- Material made-up of Cotton.

Packing	1 pc / pack
Size (Diameter x Thickness)	90 x 10 mm

MANDRELS #303

Dental Use Mandrels



Packing	12 pcs / pack
----------------	---------------

MP POWDER

Dental Polishing Material



Uses:

- Composite resin and Palladium alloys

Features:

- Used with MP BUFF, covers a whole range of polishing tasks from modifying to burnishing.
- Dust free.
- No polishing material debris deposited on the tooth neck.

Packing	1kg / pack
	1kg x 3 / box
	7kg / can

MP BUFF

Dental Polisher Buff



Uses:

- Recommended for use with MP POWDER after trimming but before final polishing.
- Hybrid Resins and Metals (using MP POWDER)

Features:

- No scattering of buff material debris.
- Removes all remaining powder clean.

Packing	1 pc / pack
Size	φ 90 x 7mm

CREAMY SAND

Dental Polishing Sand



Features:

- Sand forms like a cream making work easier and trouble-free application and polishing
- Outstanding polishing performance with brilliant luster finish.
- Cuts down polishing work time by 50%.

Packing	3kg x 2 / pack
----------------	----------------

SULFONE SAND

Dental Polishing Sand for Sulfone and Acrylic Resin



Features:

- Optimal polisher for sulfone dentures with excellent luster result.
- Outstanding polishing performance with exceptional gloss finish.
- Cuts down polishing work time by 50%.

Packing	1.5kg x 2 / pack
----------------	------------------

GLASS BEADS

Blaster Use Beads



- Uses:** #705 For Sand Blaster Use
#733 For Pencil Blaster Use

Packing	2kg / pack
Types	#705 (mesh size 149 - 250 μ)
	#733 (mesh size 44 - 88 μ)

ALUMINOUS

Blaster Use Alumina



Uniformly selected Aluminum oxide beads size for superior blasting application.

Packing	2kg / pack
Mesh Size	44 - 74 μ

POLISHING POWDER

Dental Medium Polishing Powder



Uses:
• Medium polish for metals, resin and porcelain materials.

Features:
• Substitute for polishing sand material.
• Cuts down work time and polish efficiently.
• Can be easily cleaned after use.

Packing	2kg / pack
----------------	------------

GRAZE POWDER

Dental Finish Polishing Material



Uses:
• Final polish for metal and resin materials.

Features:
• Polishing material that does not need rouge.
• Composed of fine ceramics which do not dirt hands and no effect on human body.
• All glossy polishing made easier by dissolving in water.

Packing	1.5kg / pack
----------------	--------------

SILKY SHINE

Dental Use Polisher



Uses:
• Polyamide, Polyester and Soft-Thermoplastic Resins

Features:
• Very satisfactory polishing performance is achieved when used with COTTON BRUSH.
• Liquid-type glossy finish for soft-thermoplastic resins.

Packing	30g / bottle
----------------	--------------

BLUE SHINE

Dental Final Polishing Paste



Uses: Composite Resin, Metal Alloys and Acrylic materials

Features:
• Exceptional polishing power, effortless luster and smooth finish
• Odorless results in comfortable polishing experience.
• Efficient cleaning saves polishing time

Note:
• Perform medium polishing appropriately before using BLUE SHINE.
• Too much use of polishing paste reduces polishing efficiency.

The product is water-based paste material. Water evaporates through time. Re-fill with clean water according to desired viscosity.

Packing	50g / pack
	300g / pack
	15kg / can

TIGER MULTI

Dental Medium Multi-Purpose Polishing Material



Uses: Titanium Alloy, Pure Titanium, Cobalt-Chromium Alloy, Hard-Soft Metals and Resin Polisher

Features:
• Made up of ultra fine Aluminum oxide powder that intensifies burnishing and sharpens polishing ability.
• Specially processed polisher that allows for thorough cleaning without leaving oily residue on appliances.

Packing and Size	400g / pack, 150 x 45 x 40mm
-------------------------	------------------------------

TIGER MULTI MINI

Dental Medium Multi-Purpose Polishing Material




Uses: Titanium Alloys, Pure Titanium, Cobalt-Chromium Alloy, Hard-Soft Metals and Resin Polisher

Features:
• Compact, easy to handle Tiger Multi Mini-type.
• Possible to work without touching the material directly, uses dirt-resistant plastic container.

Packing	120g
----------------	------

TIGER MULTI GOLD

Dental Medium Multi-Purpose Polishing Material



Uses: Gold Alloy, Silver Alloy and Gold-Silver-Palladium Alloy Polisher

Features:
• One product applicable both for scratch polishing and glazing of precious metal.
• Highly efficient in polishing precious metal for shiny results.
• In precious metal polishing, excessive thinning during polishing is reduced.

Packing	120g
----------------	------

ARTE SHINE

Dental Final Polishing Paste




Package	25g
----------------	-----

Type	Fine(RED) · Extra Fine(BLUE)
-------------	------------------------------


Uses: Final polishing

Range of Use: Zirconia / Porcelain/Glass ceramic /CAD/CAM /Hybrid Resin/ Hard Resin/


Image of Proper usage



【Usage Example】Steps for polishing with ARTESANO



Laboratory Equipments



LAB SCOPE S	60
LABO SCOPE S ACCESSORIES	60

LAB SCOPE S

Microscope for Dental Lab Technician



Packing	Box
Lab Scope S (body)	1 Unit
Eye Lenses (10X Magnification)	2 pcs
Mini Circle Light Joint Adapter	1 pc

Specifications

Magnification	10 X
Eye Lens	WF 10 x View 20 mm Real View 25 mm
Working Distances	120 mm
Mirror Body Formation	Straight type, rotates 360°
Mirror Body Function	Right side visibility adjustment ±5D
Eye Width Adjustment	(55 – 75) mm
Focus Adjustment	Adjustable with the flexible arm
Use Direction	Possible to fix in optional direction
Base diameter	148 mm
Relative Maximum Working Height	400 mm
Flexible Arm Length	190 mm

Usage

- Inspection of impression and plaster model surfaces
- Confirmation of margins after waxing and casting
- Examination of internal metal after casting
- Inspection of interiors and exteriors of metal bonded porcelain crowns
- Confirmation of the shifting areas on resin and porcelain
- Final inspection of finished prosthesis

Features

- Compact size and lightweight, easy to handle and requires little bench space
- Flexible neck allows angle adjustment, direction and height
- Protective Lens Cover supplied

LABO SCOPE S ACCESSORIES



LED CIRCLE LIGHT (100V)

Packing	
Bulb Light	1 pc
Bulb Light Holder	1 pc



MINI CIRCLE LIGHT JOINT ADAPTER

Packing	
Joint Adapter	1 pc

*Customers without the joint adapter need to purchase one order to attach Mini Circle Light.



LENSES

Packing	
WF5 Lenses (Magnification Power: 5X)	2 pc / set
WF20 Lenses (Magnification Power: 20X)	2 pc / set

Attachment



What is KUGEL HOOK?	62
KUGEL HOOK	62

Attachment

What is KUGEL HOOK?

KUGEL HOOK offers solution for two aspects of dental application, namely tooth lose and denture mechanical stability. There are situations where decaying or severely damaged tooth has turned beyond repair. In this instance, tooth is extracted and a denture is consequently replaced. On the denture part, a variety of products are available which promotes mechanical stability of the denture relative to its surrounding mouth and gum. These include abutments, clasps and braces. This is the conventional process of resolving issues from tooth lose to denture replacement.


On the other hand, **KUGEL HOOK** has been conceptualized in order to alternatively abridge tooth lose and denture mechanical stability relationship. Along the process, instead of losing the tooth – **KUGEL HOOK** invokes utilization of its base and thence transforming into a denture hook. In this way, without losing the tooth completely, the gum integrity and natural teeth alignment are preserved.



KUGEL HOOK Attachment

Packing	
Male	1pc / box, 5pcs / pack
Male Holder	1pc / pack
Female	1pc / box (Setting Film x 2) 5pcs / pack (Setting Film x 10)

KUGEL HOOK is composed of male and female parts. The male part is used as the bolt impression of the tooth base for metal casting. The metal casting is cemented into the excavated tooth base. The plastic female part is precisely affixed in the interior part of the denture using self-curing resin, as a socket, where the bolt is to be attached. It acts as a bolt-and-socket device between supposedly gum and denture and therefore guarantees denture mechanical stability against grinding and chewing.



KUGEL HOOK portrays a semi-implant conservative approach addressing the matter over denture mechanical stability without sacrificing the tooth of concern entirely for a much more economical and faster recovery than any conventional implant technique.



PRODUCT CATALOG



YAMAHACHI DENTAL

54-1, Ochigara, Nishiura-Cho
Gamagori-City, Aichi - JAPAN

TEL : +81-533-57-7121
FAX : +81-533-57-1764

box@yamahachi-dental.co.jp
www.yamahachi-dental.co.jp/en/