

# GLOSSY COAT COLOR

Light Curing Naturalizing Glossing Agent



Can be used in the lab (Extraoral) or chair side (Intraoral).  
Quick and full control of the shade of resin based prosthesis.



YAMAHACHI DENTAL

# GLOSSY COAT COLOR

## Color Resistance


Excellent color resistance after curing.



Immersion in 2% coffee solution at 37°C for 1 week

## Abrasion Resistance


Superior abrasion resistance after curing.



Toothbrush Abrasion Test after 100,000 cycles

## Adhesion


Excellent adhesion to many types of resin based materials.




Flexural Bond Strength

## Directions of use -Chair side (Intraoral) -

**\*This product is for resin based prosthesis. It cannot be used on natural teeth.**



## Directions of use -Lab side (Extraoral) -



## Example of Use

After Surface Roughing



After using Clear liquid only



Brown: Clear = 1:1 after use



After using Brown liquid only



After characterisation, the surface is coated with clear liquid.

## Product Information



**Packing** Clear 8ml/pc  
Brown 5ml/pc  
Wavelength 365-410 nm

Being distributed by:



**YAMAHACHI DENTAL**

54-1, OCHIGARA, NISHIURA-CHO, GAMAGORI, AICHI-PREF., 443-0105, JAPAN

TEL : 81-533-57-7121 FAX : 81-533-57-1764

Website <https://yamahachi-dental.co.jp/en/>

E-mail [box@yamahachi-dental.co.jp](mailto:box@yamahachi-dental.co.jp)