

TrinDy

Cast

3D Printing Resin

TPO Free



Packing : 1,000 g bottle / 500 g bottle

Color :



Yellow

Contact us

Designed for dental cast fabrication. Excellent printing accuracy and clean burn out characteristics on micro-wax are taken into account in the resin structural design, while not only meeting to the basic requirement of performance as fabricated dental casts.

- Patterns for casting and pressing
- Crown and bridge
- Removable partial denture frameworks



| Property | Procedure | TrinDy Cast |
|-------------------|------------|-------------|
| Flexural Strength | MPa ISO178 | 20 |
| Flexural Modulus | MPa ISO178 | 640 |
| Shore Hardness | - ISO48-4 | ≥D65 |
| Viscosity | mPa·s | 300-400 |

Printer : ASIGA MAX (ASIGA)

Second curing : Machine(Otoflash)/Condition(1 × 2000 flashes)

Items mentioned in this information are for customers reference and limited to certain examples although they are based on the most reliable data obtained by our company. It is hoped that customers make the best use of these products after careful examination.