




YAMAHACHI DENTAL

# TrinDy®

## Mask

3D Printing Resin

TPO Free



Packing : 1,000 g bottle

Color :



Pink

Contact us

Designed for gingiva mask fabrication. Its excellent printing accuracy, elasticity and sufficient tear strength enable operators to look and feel the gum tissue on simulation.

- Gingiva masks
- Soft tissue for implant models



Pink

Property		Procedure	TrinDy Mask
Tensile Strength	MPa	ASTM D412	2
Elongation at Break	%	ASTM D412	110
Shore Hardness	-	ISO48-4	≥A65
Viscosity	mPa·s		850-950

Printer : ASIGA MAX (ASIGA)

Second curing : Machine(Otoflash)/Condition(1 × 2000 flashes)

Items mentioned in this information are for customers reference and limited to certain examples although they are based on the most reliable data obtained by our company. It is hoped that customers make the best use of these products after careful examination.

