



YAMAHACHI DENTAL MFG., CO.

MANUFACTURER, IMPORTER & EXPORTER OF DENTAL PRODUCTS

54-1 Ochigara, Nishiura-Cho, Gamagori City, Aichi Prefecture, 443-0105 Japan

Tel: +81-533-57-7121 Fax: +81-533-57-1764

URL: <http://www.yamahachi-dental.co.jp> e-mail: box@yamahachi-dental.co.jp

December 27, 2024

Notification of change in tooth plates design

To Whom It May Concern:

This is to inform you all just in case that the portion of the plastic plates for posterior teeth where the artificial teeth are placed has been changed from Tape to Wax since the May 2024 production.

This is to prevent the artificial tooth from falling out of the plate, and the following design changes have also occurred at the same time. And the same changes also occurred for the plates for anterior teeth.

1. Removal of QR code on surface

This QR code is for our factory management and does not affect our customers in any way.

The factory operates the UDI code on the reverse side instead.

*An example of the UDI code attached to the back of the product




2. Removal of Lot. notation on surface

Lot numbers can be found on a sticker seal on the back instead.

*Actual images

Old plates:



New plates: Can hold even large posterior teeth (size# 36) securely



Please be advised that our shipments will be mixed with the old and the new for a while.

We ask for your understanding.

Sincerely yours,

遠山昌志



Masashi Toyama / President