

COMPOSITE RESIN TEETH EFUCERA-A & EFUCERA-P MOULD CHART

EFUCERA-A


▲ T4	40.2 7.7 10.3	■ SS6	46.0 8.9 10.6
▲ T5	44.4 8.4 11.0	■ SS7	48.4 9.2 11.0
▲ T6	46.2 8.9 11.5	■ C4	42.2 7.8 10.7
▲ T7	49.0 9.3 11.9	■ C5	45.0 8.3 11.3
■ S4	43.4 7.9 10.4	■ C6	47.2 8.7 11.6
■ S5	45.4 8.3 10.9	■ C7	49.4 9.1 12.2
■ S6	48.4 8.8 11.4	■ LA4	4.9 33.2 9.5
■ S7	50.8 9.2 12.0	■ LA5	5.2 34.6 10.6
■ SS4	40.6 7.8 9.5	■ LA6	5.5 37.4 10.9
■ SS5	43.6 8.3 10.1	■ LA7	5.9 39.4 11.0



EFUCERA-A



EFUCERA-P



COMBINATION TABLE

ANTERIOR		POSTERIOR
EFUCERA-A		EFUCERA-P
UPPER	LOWER	
T4	LB4	28
T5	LB5	30
T6	LB6	32
T7	LB7	34(36)
S4	LA4	30
S5	LB6	30
S6	LB7	32
S7	LA7	34(36)
SS4	LB4	28
SS5	LB5	30
SS6	LA6	32
SS7	LB7	34(36)
C4	LA4	30
C5	LA5	32
C6	LA6	32
C7	LA7	34(36)

Registered in the excad and 3 shape libraries

*Cusp Angle 20 degrees

SHADES A1 ,A2 ,A3 ,A3.5 ,A4

DISTRIBUTOR



YAMAHACHI DENTAL

54-1, OCHIGARA, NISHIURA-CHO, GAMAGORI, AICHI-PREF., 443-0105, JAPAN

TEL : 81-533-57-7121 FAX : 81-533-57-1764

Website <https://yamahachi-dental.co.jp/en/>

E-mail box@yamahachi-dental.co.jp